MODELCAST S

EN ISO 15912
03-2014

NAVODILA ZA UPORABO
WORKING INSTRUCTIONS
GEBRAUCHSANWEISUNG
NÁVOD NA ZPRACOVÁNÍ
UPUTA ZA KORIŠTENJE
ИНСТРУКЦИИ ПО ПРИМЕНЕНИЮ
UPUTSTVO ZA UPOTREBU
UPUTSTVO ZA KORIŠĆENJE
УПАТСТВО ЗА КОРИСТЕЊЕ

Proizvajalec/Producer/Hersteller
SI - 3000 CELJE · Opekarniška cesta 26
Tel.: +386 (0)3 425-62-00 · Fax: +386 (0)3 425-62-10
E-mail:info@interdent.cc · http://www.interdent.cc

CZ
INTERDENT s.r.o.
Foerstrova 12, Strašnice
CZ-10000 Praga
tel. +420/274 814 -404
fax +420/274 820 130

HR
INTERDENT d.o.o
Vinogradski odvojak 2d
HR-10431 Sveta Nedelja
tel. +385/1 3873 644
fax +385/1 3873 617
mail: interdent@interdent.hr

SRB
INTERDENT d.o.o.
Nehruova 44
SRB-11070 Novi Beograd
tel. +381/11 217 53 74
fax +381/11 3 01 75 33

MK
INTERDENT d.o.o.
Bojmija 1-12M
MK-1000 Skopje
tel. +389/224 628 55
fax +389/224 634 12

SK
INTERDENT s.r.o.
Liščie údolie 57
SK-84231 Bratislava
tel. +421/254 413 697
fax +421/254 413 626

NAVODILA ZA UPORABO MODELCAST S — SI

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WORKING INSTRUCTIONS MODELCAST S — GB

Modelcast S is extremely fine grained, temperature resistant phosphate investment material used for making partial dentures with metal base for speed or conventional heating method. It can be duplicated only in silicone. For good mixture, good chemical binding and smooth cast the investment should be mixed 15 seconds by hand and 60 seconds by vacuum.

- Before you mix the powder and liquid the mixing bowl has to be clean and dry.
- When you use the silicon for duplicating and you place it under the pressure, you must continue the complete process under pressure.
- Use different mixing bowl for investment and stone. Take care that the investment does not come into the contact with stone.
- We recommend using degreasing agent for wax parts. (Interwaxit). Interwaxit is necessary well-dry after use.
- Keep furnace door closed for 10 minutes after insertion of the mould!

WORKING INSTRUCTIONS:
- First pour liquid and then add powder.
- Manual mixing time 15 s
- Vacuum mixing time 60 s
- Duplicate model in silicone.
- Prepare the duplicated model (The base of the model should have a minimum thickness of 1 cm) –
- Remove the model from silicon after 30-40 minutes.
- Dip/Spray the duplicated model: MODELHART, MODELHART SPRAY
- Do the wax modelation
- Invest into the mould (see the mixing ratio).
- Scratch the surface of the mould before place it in the furnace

MIXING RATIO POWDER/ LIQUID (DUPLICATION IN SILICONE)
200g / 40ml.....model
400 g/ 90ml.....mould

LIQUID MIX (MIXING RATIO):

SILIKON		
	Expasol	DISTILLED WATER
PARTIAL DENTURE	Model: (75-80%) Mould: (50-70%)	Model: (25-20%) Mould: (50-30%)
COMBI TECHNIQUE	Model: (80-85%) Mould: (50-70%)	Model:(20-15%) Mould: (50-30%)

MIXING RATIO:

	Expasol	Distilled Water	Indication
Model 200g/40ml	75% (30ml)	25% (10ml)	Partial denture
	80% (32 ml)	20% (8 ml)	Partial denture, combi work
	85% (34 ml)	15% (6 ml)	combi work
Mould 400g/90ml	60% (54 ml)	40% (36 ml)	Partial denture, combi work

SPEED HEATING METHOD
Place the mould in the furnace at 900- 950°C 17-20 minutes after first contact of liquid and powder.
Holding time minimum .60 minutes – (depending on size and number of moulds).

CONVENTIONAL HEATING METHOD:

Temperature	Rate of heating	Holding time
Room temperature– 270°C	5°C/min	40min
270°C-580°C	7°C/min	30min
580°C – 950°C	9°C/min	60 min

Holding time minimum .60 minutes – (depending on size and number of moulds).
Casting (I-MG, I-MG FH, I-MG EKO)
Let the mould cool down to the room temperature, do not force the cooling.
After cooling to the room temperature take the cast out of the investment.

WARRANTY!
Whether given verbally, in writing or by practical instructions, our recommendations for use are based upon our own experience and trials and can only be considered as standard values.
Our products are subject to further development. Therefore changes in construction and composition are reserved.

STORAGE: at room temperature
INVESTING TEMPERATURE: 19°C – 23°C
PERSONAL PROTECTION: We recommend wearing half facemask with filter FFP1.

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