# INTERVEST K + B SPEED

**EN ISO 15912**
05-2013

Navodila za uporabo / Working instructions / Gebrauchsanleitung /
Instrucciones de uso / Návod k použití / Οδηγιες χρησης / Upute za uporabu
Дополнительные инструкции по применению / Uputstvo za upotrebu
Instrukcije za upotrebu / Упатство за употреба

**Proizvajalec/Producer/Hersteller**
SI - 3000 CELJE · Opekarniška cesta 26
Tel.: +386 (0)3 425-62-00 · Fax: +386 (0)3 425-62-10
E-mail:info@interdent.cc · http://www.interdent.cc

**CZ**
INTERDENT s.r.o.
Foerstrova 12, Strašnice
CZ-10000 Praga
tel. +420/274 814 -404
fax +420/274 820 130

**HR**
INTERDENT d.o.o
Vinogradski odvojak 2d
HR-10431 Sveta Nedelja
tel. +385/1 3873 644
fax +385/1 3873 617
mail: interdent@interdent.hr

**SRB**
INTERDENT d.o.o.
Nehruova 44
SRB-11070 Novi Beograd
tel. +381/11 217 53 74
fax +381/11 3 01 75 33

**MK**
INTERDENT d.o.o.
Bojmija 1-12M
MK-1000 Skopje
tel. +389/224 628 55
fax +389/224 634 12

**SK**
INTERDENT s.r.o.
Líšće údolie 57
SK-84231 Bratislava
tel. +421/254 413 697
fax +421/254 413 626

## NAVODILA ZA UPORABO INTERVEST K + B SPEED

SI

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SKLADIŠČENJE: pri sobni temperaturi
TEMPERATURA VLAGANJA: 19°C - 23°C

## WORKING INSTRUCTIONS INTERVEST K + B SPEED

GB

**Shock heat or conventionally heatable, phosphate-bonded investment material.**

Mixing ratio powder/liquid : 160g (1 bag)/40 ml liquid

**HIGHER EXPASOL CONCENTRATION → MORE EXPANSION** → BIGGER CROWN → If the crown is to big, the concentration of Expasol should be reduced.

1. Wax-up using acrylics + 5% EXPASOL
2. Mixing bigger quantities of investment - 5% EXPASOL
3. Using of metal ring 6X or 9X is recommended. Always use two liners when using ring number 6 and 3 liners when using ring number 9 to provide good expansion.
4. Do not wet ceramic liner with water.
5. For good mixture, good chemical binding and smooth cast the investment should be mixed **15 seconds by hand and 60 seconds by vacuum**.
6. Use different mixing bowl for investment and stone. Take care that investment does not come into contact with stone.

**RAPID FIRING:** important notes
1. Roughen the mould surface and place the mould with the funnel former facing down and without direct contact to the bottom or the walls in the furnace. (Use distance holder or ceramic plate.)
2. Place the moulds in the furnace, which has been heated up to final temperature, 15-20 minutes after first contact of liquid and investment. The final temperature depends on the type of alloy ( eg. I bond 02 and I bond NF 900 $^0$C ).
3. Maintain the final temperature at least 30 minutes for mould number 1, 40 minutes for the mould number 3, 60 minutes for the mould number 6 and 90 minutes for the mould number 9.
4. Keep furnace door **10 minutes closed** after insertion of the mould!

**CONVENTIONAL BURNOUT:** If you do not place the mould in the furnace after 15 -20 minutes you should use conventional burnout.

**Furnaces with conventional control**
Place mould in cold furnace and heat it up to 270°C. Maintain at this temperature 45-60 minutes. Heat up to final temperature and maintain it 30 - 90 minutes, depending on the size of the mould.

**Furnaces with computer control**

| Temperature | Rate of heating | Holding time |
|---|---|---|
| Room temperature - 270°C | 7°C/min | 40min |
| 270°C-580°C | 7°C/min | 30min |
| 580°C - final temperature (I-Bond 02, I-Bond nf 900$^0$C ) | 9°C/min | 30min |

**Warranty!** Whether given verbally, in writing or by practical instructions, our recommendations for use are based upon our own experience and trials and can only be considered a standard values. Our product are subject to further development Therefor changes in construction and composition are reserved.

**Personal protection:** We recommend wearing half facemask with filter A1.

STORAGE: at room temperature
INVESTING TEMPERATURE: 19°C - 23°C

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